



# RESTORIGIN™

## Human Amniotic Membrane (HAM)

Its Effective Use in the Management of Chronic or  
Acute Diabetic Wounds and Venous Leg Ulcers - A Case Study

Milla Zakharova, PhD, Brent Hall, Mark Schallenberger, Keith Bangart, DPM,  
Dan Bangart, DPM, Steven Moore, DPM, Jeffrey Thomas, DPM

# Case study report of chronic non-healing foot ulcers treated with dehydrated human amniotic membrane sheet

## Background

Chronic non-healing wounds should initially be treated in accordance with standard care protocols<sup>1,2</sup>, however, the wounds that did not achieve 53% area reduction by 4 weeks healed in only 9% of the cases ( $p < 0.001$ ).<sup>3</sup> Current treatment options are often insufficient to achieve complete wound closure. Dehydrated human amniotic membrane (HAM) is non-immunogenic biological material with anti-microbial, anti-inflammatory and anti-fibrotic properties that function as a scaffold to facilitate tissue regeneration.<sup>4</sup>

## Study Objective

The purpose of this study was to evaluate the efficacy of a sterile, dehydrated HAM sheet (Restorigin™) for treating non-healing wounds of the lower extremities.

## Study Design

### Patient population and inclusion criteria

- 10 patients with Ulcer Wagner Grades 1 and 2:
  - Five patients with Type I or Type II Diabetes mellitus;
  - One patient with peripheral vascular disease
- Failure to improve after 4 weeks standard of care

### Treatment

- Restorigin™ weekly or bi-weekly application

## Results

- The baseline ulcer area was  $2.3 \pm 0.9 \text{ cm}^2$
- Median number of Restorigin™ applications was 3.5 (from 1 to 7)
- In two patients (20%), wounds were completely healed only two weeks after initial application of the Restorigin™ membrane sheet
- On average, 94.6% of wound area reduction was observed after 8 weeks of HAM therapy.




Figure 1. Individual wound closure results for each case.



Figure 2. Healing progression is shown for a 42-year-old female patient from before the treatment (A) to complete closure after 4.5 weeks of HAM therapy and a single graft application (C).




Figure 3. An average of 94.6% of wound area reduction was observed after 8 weeks of HAM sheet therapy.

## Conclusion

The study demonstrated that human amniotic membrane (HAM) graft therapy is a safe and effective treatment for chronic non-healing ulcers.

1. Frykberg R, et al. J Foot Ankle Surg. 2006 Sep-Oct;45(5 Suppl): S1-66.
2. Bakker K, et al; Diabetes Metab Res Rev. 2012 Feb;28 Suppl 1:225-31.
3. Sheehan P, Jones P, et al. Diabetes Care. 2003 Jun; 26(6):1879-82.
4. Hao Y, Ma D, et al. Cornea. 2000 May;19(3):348-52.